



Healthcare Range












www.reposefurniture.co.uk

Seat
Cushions



Seat Cushions

Product not available on small chairs

			Multi C-Air	Multi-Flex	Meirose	Flexi Porter	C-air	Brooklyn	Madison	Multi Bari	Bariatric	Stargazer	Harlem	Boston	Stamford	Washington	Nashville	Arden	Haven
 <p>Easiform</p> <p>For standard seating, our softest foam cushion with high quality Reflex foam filler for gentle comfort.</p>			✓		✓		✓			✓	✓	✓			✓	✓		✓	
 <p>Firmaform</p> <p>For our standard seating, a cushion that provides a more supportive and firm seat with a high-quality filler from Reflex, the foam specialists.</p>			✓		✓		✓			✓	✓							✓	
 <p>Foamaform</p> <p>For those who have good core stability and posture, are sitting for more than one hour at a time and require low to medium pressure relief combined with a high level of comfort. This cushion has different layers of high-density foam, offering a supportive seat base which allows the user's bottom to be immersed into the cushion. It provides stability, a large surface area reducing pressure levels and a high level of personal comfort.</p>														✓					
 <p>Viscoform</p> <p>For those who have good posture and mobility but are seated for greater than one hour at a time and require pressure relief and supported comfort. This cushion has a high-density foam base which offers a supportive and firm seat base with a top layer of memory foam to allow immersion into the cushion, giving stability and a large surface area to reduce pressure levels.</p>			✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓
 <p>Profile Viscoform</p> <p>The profile version of the Viscoform cushion provides aommel, which can support those with weak trunk or lower limb muscles who are prone to sliding down and forward in their chair or for those with lower limb spasticity to help maintain hips in the abducted position to create a good seating position.</p>			✓		✓		✓	✓											
 <p>Liquiform</p> <p>For those who have good posture and mobility but sit for long periods of time and require pressure relief and supported comfort. Recommended for people with limited muscle bulk or at risk of excessive weight loss. This cushion has a high-density foam which offers a supportive and firm seat base. The foam is topped with a fluid liquid gel in a soft silicone cover which provides a high level of pressure relief. The liquid gel helps to dissipate heat and constantly adapts to the user's movement.</p>			✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓		
 <p>Coolform</p> <p>For those sitting for long periods with specific comfort issues such as a high risk of pressure sores. This cushion uses the latest medical grade silicone gel technology allowing pressure to be evenly distributed. The gel conforms to the shape of the body and equalises to body temperature. The gel cells allow air to flow easily and Coolform is totally odourless and hypoallergenic.</p>			✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓		✓	✓	
 <p>Profile coolform</p> <p>The profile version of the Coolform cushion provides aommel, which can support those with weak trunk or lower limb muscles who are prone to sliding down and forward in their chair or for those with lower limb spasticity to help maintain hips in the abducted position to create a good seating position.</p>			✓		✓		✓	✓											
 <p>Celliform</p> <p>For those sitting for long periods with specific comfort issues such as a high risk of pressure sores. Gelacell has a unique cell matrix that distributes pressure across the contact area, as it contours and yields in response to weight redistribution. Gelacell absorbs body heat, reducing the risk of creating an environment that could increase pressure risks, the grid design also encourages air circulation. Gelacell is clinically tested and passes all dermatological tests for skin contact.</p>			✓	✓	✓	✓	✓	✓					✓	✓			✓		
 <p>Airform</p> <p>For those who currently have or who are at high risk of pressure sores. Users need to have good pelvic stability and core trunk strength to accommodate the air pressure changes. The Airform has three manual air cushions under a layer of memory foam. The cushion provides all the benefit of an air cushion without the sensation of alternating changes in air pressure.</p>			✓		✓		✓	✓						✓					✓
 <p>Dynaform</p> <p>For those at very high risk of developing pressure sores or a history of sores where the highest priority is to prevent further pressure damage. The Dynaform cushion is a dynamic alternating air cushion, which constantly changes the air levels within its air cells to recreate the effect of normal body movements, which provides constant pressure redistribution to prevent pressure build up.</p>			✓		✓		✓	✓	✓								✓		



www.reposefurniture.co.uk



Seat
Cushions